


CD73 Antibody [clone RM431] (R20446)

Catalog No.	Formulation	Size
R20446-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM431
Purity	Protein A purified from animal origin-free supernatant
UniProt	P21589
Localization	Cell surface, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:100 -1:200 Western Blot : 1:500-1:1000
Limitations	This recombinant CD73 antibody is available for research use only.



Western blot testing of human MDA-MB-231 cell lysate with recombinant CD73 antibody at 1:1000. Predicted molecular weight: 60-70 kDa.



IHC staining of FFPE human liver cancer tissue with recombinant CD73 antibody at 1:200.

Description

This antibody reacts to human CD73 (also called ecto-5'-nucleotidase, E5NT, NT, NT5, NTE, eN and eNT.)

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant CD73 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A recombinant human CD73 protein was used as the immunogen for the CD73 antibody.

Storage

Store the recombinant CD73 antibody at -20oC.