


CD8a Antibody [clone RM397] (R20413)

Catalog No.	Formulation	Size
R20413-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM397
Purity	Protein A purified from animal origin-free supernatant
UniProt	P01732
Localization	Cell surface
Applications	Immunohistochemistry (FFPE) : 1:100-1:200 Western Blot : 1:100-1:200
Limitations	This recombinant CD8a antibody is available for research use only.



Western blot testing of human T-cell lysate with recombinant CD8a antibody at 1:100 dilution. Expected molecular weight: 26-40 kDa depending on glycosylation level.



IHC staining of FFPE human appendix tissue with recombinant CD8a antibody at 1:100.

Description

This antibody reacts to the alpha chain of human CD8.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant CD8a antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to human CD8a was used as the immunogen for the recombinant CD8a antibody.

Storage

Store the recombinant CD8a antibody at -20oC.