

CD90 Antibody / Thy1 [clone AF-9] (V2892)


Catalog No.	Formulation	Size
V2892-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V2892-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V2892SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug



Citations (3)

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	AF-9
Purity	Protein G affinity chromatography
UniProt	P04216
Localization	Cell surface
Applications	Flow Cytometry : 0.5-1ug/million cells Immunofluorescence : 0.5-1ug/ml
Limitations	This CD90 antibody is available for research use only.



SDS-PAGE Analysis of Purified, BSA-Free CD90 Antibody (clone AF-9). Confirmation of Integrity and Purity of the Antibody.

