

CDH10 Antibody / Cadherin 10 [clone 28C85] (RQ8643)

Catalog No.	Formulation	Size
RQ8643	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ug

Recombinant RABBIT MONOCLONAL

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28C85
Purity	Affinity chromatography
UniProt	Q9Y6N8
Applications	Western Blot : 1:500-1:2000
Limitations	This CDH10 antibody is available for research use only.



Western blot testing of human U-87 MG cell lysate with CDH10 antibody. Predicted molecular weight ~88 kDa but may be observed at higher molecular weights due to glycosylation.

Description

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. [UniProt]

Application Notes

Optimal dilution of the CDH10 antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human Cadherin 10 protein was used as the immunogen for the CDH10 antibody.

Storage

Store the CDH10 antibody at -20oC.