

CDK9 Antibody (C-Terminal Region) (R36496)

Catalog No.	Formulation	Size
R36496-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Dog, Human, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q99J95
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This CDK9 antibody is available for research use only.



Western blot testing of mouse NIH3T3 cell lysate with CDK9 antibody at 2ug/ml.
Predicted molecular weight: ~43 kDa.

Description

Alternate name(s) for this target: Cyclin-dependent kinase 9; C-2k; CDC2L4; CTK1; PITALR; TAK; CDC2-related kinase, cell division cycle 2-like protein kinase 4; cell division protein kinase 9; serine/threonine protein kinase PITALRE; tat-associated kinase complex catalytic subunit

Application Notes

Optimal dilution of the CDK9 antibody should be determined by the researcher.

Immunogen

Amino acids KGSQITQQSTNQSRN were used as the immunogen for this CDK9 antibody.

Storage

Aliquot and store the CDK9 antibody at -20oC.