

CDKN2A Antibody (pancreatic isoform specific) (R34337)

Catalog No.	Formulation	Size
R34337-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1029
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CDKN2A antibody is available for research use only.



Western blot testing of HeLa lysate with CDKN2A antibody at 1ug/ml. Predicted molecular ~12kDa, observed here at ~17kDa.

Description

Additional name(s) for this target protein: Cyclin-dependent kinase inhibitor 2A; p16INK4; p16

Application Notes

Optimal dilution of the CDKN2A antibody should be determined by the researcher.

1. This CDKN2A antibody detects the pancreatic form of the protein.

Immunogen

Amino acids ILRKKGRLPGEFSE were used as the immunogen for this CDKN2A antibody.

Storage

Aliquot and store the CDKN2A antibody at -20oC.