


CDT1 Antibody (R33167)

Catalog No.	Formulation	Size
R33167-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81620
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CDT1 antibody is available for research use only.



Western blot testing of human ovary lysate with CDT1 antibody at 0.05ug/ml. Predicted molecular weight ~60 kDa, commonly observed at 60-70 kDa.

Description

Additional name(s) for this target protein: Chromatin licensing and DNA replication factor 1

Application Notes

Optimal dilution of the CDT1 antibody should be determined by the researcher.

Immunogen

Amino acids ARLAHQTRAEEGL were used as the immunogen for this CDT1 antibody.

Storage

Aliquot and store the CDT1 antibody at -20oC.