


CEA Antibody / Carcinoembryonic Antigen [clone RM326] (R20351)

Catalog No.	Formulation	Size
R20351-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM326
Purity	Protein A purified from animal origin-free supernatant
UniProt	P06731
Localization	Cytoplasmic and luminal surface
Applications	Immunohistochemistry (FFPE) : 1:1000-1:2000 Western Blot : 1:1000-1:2000 ELISA : 1:1000-1:5000
Limitations	This recombinant CEA antibody is available for research use only.



Western blot testing of human MCF7 cell lysate with recombinant CEA antibody at 1:1000. Expected molecular weight: 80–200 kDa depending on glycosylation level.



IHC staining of FFPE human colon tissue with recombinant CEA antibody at 1:2000.

Description

This antibody reacts to human CEA (Carcinoembryonic Antigen, CD66e).

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant CEA antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Native Carcinoembryonic Antigen protein purified from human liver tissue was used as the immunogen for the recombinant CEA antibody.

Storage

Store the recombinant CEA antibody at -20oC.