


Chromogranin A Antibody / CHGA [clone RM385] (R20402)

Catalog No.	Formulation	Size
R20402-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM385
Purity	Protein A purified from animal origin-free supernatant
UniProt	P10645
Localization	Finely granular cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:5000-1:10000 Western Blot : 1:1000-1:2000
Limitations	This recombinant Chromogranin A antibody is available for research use only.




Western blot testing of human brain tissue with recombinant Chromogranin A antibody at 1:1000. Predicted molecular weight ~51 kDa but may be observed at higher molecular weights due to glycosylation.



IHC staining of FFPE human pancreas tissue with recombinant Chromogranin A antibody at 1:2500.



IHC staining of FFPE human colon tissue with recombinant Chromogranin A antibody at 1:10,000.



IHC staining of FFPE human medullary thyroid carcinoma tissue with recombinant Chromogranin A antibody at 1:2500.

Description

This antibody reacts to human Chromogranin A.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Chromogranin A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the N-terminus of human Chromogranin A was used as the immunogen for the Chromogranin A antibody.

Storage

Store the recombinant Chromogranin A antibody at -20°C.