


CHUK Antibody (R33791)

Catalog No.	Formulation	Size
R33791-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1147
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CHUK antibody is available for research use only.



Western blot testing of HeLa nuclear lysate with CHUK antibody at 1ug/ml. Predicted molecular weight: ~85kDa.

Description

Additional name(s) for this target protein: Conserved helix-loop-helix ubiquitous kinase

Application Notes

Optimal dilution of the CHUK antibody should be determined by the researcher.

Immunogen

Amino acids EHDHSLSCVVTPQD were used as the immunogen for this CHUK antibody.

Storage

Aliquot and store the CHUK antibody at -20oC.