

cIAP2 Antibody (F49443)


Catalog No.	Formulation	Size
F49443-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49443-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q13489
Localization	Cytoplasmic
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This cIAP2 antibody is available for research use only.



cIAP2 antibody western blot analysis in A549 lysate



cIAP2 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human brain tissue.

Description

BIRC3 is a member of a family of proteins that inhibits apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. The encoded protein inhibits apoptosis induced by serum deprivation but does not affect apoptosis resulting from exposure to menadione, a potent inducer of free radicals. The amino acid sequence predicts three baculovirus IAP repeat domains and a ring finger domain.

Application Notes

Titration of the cIAP2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 101-131 from the human protein was used as the immunogen for this cIAP2 antibody.

Storage

Aliquot the cIAP2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.