


CIRBP Antibody (R35190)

Catalog No.	Formulation	Size
R35190-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1153
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This CIRBP antibody is available for research use only.



Western blot testing of mouse testis lysate with CIRBP antibody at 0.05ug/ml. Both the expected ~18kDa band and the additional ~23kDa band are blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Cold inducible RNA binding protein

Application Notes

Optimal dilution of the CIRBP antibody should be determined by the researcher.

Immunogen

Amino acids QAGKSSDNRSR were used as the immunogen for this CIRBP antibody.

Storage

Aliquot and store the CIRBP antibody at -20oC.