


CK20 Antibody / KRT20 / Cytokeratin 20 (R34922)

Catalog No.	Formulation	Size
R34922-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54474
Applications	Western Blot : 0.01-0.03ug/ml Enzyme ImmunoAssay : used as capture for recombinant protein ELISA (peptide) LOD : 1:128000
Limitations	This CK20 antibody is available for research use only.



Western blot testing of human colon lysate with CK20 antibody at 0.03ug/ml. Predicted molecular weight: ~49kDa.

Description

Additional name(s) for this target protein: KRT20, keratin 20, CK-20, Cytokeratin 20

Application Notes

Optimal dilution of the CK20 antibody should be determined by the researcher.

Immunogen

Amino acids NSRHTVNYGSDLT were used as the immunogen for this CK20 antibody.

Storage

Aliquot and store the CK20 antibody at -20oC.