


Claudin 2 Antibody (F49581)

Catalog No.	Formulation	Size
F49581-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49581-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P57739
Applications	Western Blot : 1:1000
Limitations	This Claudin 2 antibody is available for research use only.



Description

Members of the claudin protein family, such as CLDN2, are expressed in an organ-specific manner and regulate the tissue-specific physiologic properties of tight junctions.

Application Notes

Titration of the Claudin 2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 190-220 from the human protein was used as the immunogen for this Claudin 2 antibody.

Storage

Aliquot the Claudin 2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.