


COG7 Antibody (R34676)


Catalog No.	Formulation	Size
R34676-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	91949
Localization	Cytoplasmic
Applications	Western Blot : 0.03-0.1ug/ml IHC (FFPE) 5ug/ml : ELISA (peptide) LOD : 1:32000
Limitations	This COG7 antibody is available for research use only.



Western blot testing of MOLT4 lysate with COG7 antibody at 0.03ug/ml. Predicted molecular weight: ~86 kDa.



IHC testing of FFPE human testis tissue with COG7 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Component of oligomeric golgi complex 7

Application Notes

Optimal dilution of the COG7 antibody should be determined by the researcher.

Immunogen

Amino acids KTRPEDYRQVSK were used as the immunogen for this COG7 antibody.

Storage

Aliquot and store the COG7 antibody at -20oC.