


COL4A3BP Antibody (R34612)

Catalog No.	Formulation	Size
R34612-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10087
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This COL4A3BP antibody is available for research use only.



Western blot testing of human heart lysate with COL4A3BP antibody at 1ug/ml.
Predicted molecular weight: ~70kDa.

Description

Additional name(s) for this target protein: Collagen, type IV, alpha 3 (Goodpasture antigen) binding protein

Application Notes

Optimal dilution of the COL4A3BP antibody should be determined by the researcher.

Immunogen

Amino acids KREDWQKRLDKETEK were used as the immunogen for this COL4A3BP antibody.

Storage

Aliquot and store the COL4A3BP antibody at -20oC.