


Collagenase 3 Antibody / MMP13 (RQ6656)

Catalog No.	Formulation	Size
RQ6656	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P45452
Applications	Western Blot : 1-2ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This Collagenase 3 antibody is available for research use only.



Western blot testing of human HepG2 cell lysate with Collagenase 3 antibody. Predicted molecular weight ~54 kDa.

Description

Collagenase 3 is an enzyme that in humans is encoded by the MMP13 gene. This gene encodes a member of the peptidase M10 family of matrix metalloproteinases (MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature protease. This protease cleaves type II collagen more efficiently than types I and III. It may be involved in articular cartilage turnover and cartilage pathophysiology associated with osteoarthritis. Mutations in this gene are associated with metaphyseal anadysplasia. This gene is part of a cluster of MMP

genes on chromosome 11.

Application Notes

Optimal dilution of the Collagenase 3 antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids K249-C471) was used as the immunogen for the Collagenase 3 antibody.

Storage

After reconstitution, the Collagenase 3 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.