


Contactin 1 Antibody / GP135 (R33893)


Catalog No.	Formulation	Size
R33893-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1272
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Contactin 1 antibody is available for research use only.



Western blot testing of human cerebellum [A], mouse brain [B] and rat brain [C] lysates with Contactin 1 antibody at 0.3ug/ml. Predicted molecular weight: ~113/135-150kDa (unmodified/glycosylated).



Western blot of pig brain lysate with Contactin 1 antibody at 0.3ug/ml. Predicted molecular weight: ~113/135-150kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: CNTN1; neural cell surface protein F3

Application Notes

Optimal dilution of the Contactin 1 antibody should be determined by the researcher.

Immunogen

Amino acids ENIRGKDKHQAR (internal region) were used as the immunogen for this Contactin 1 antibody.

Storage

Aliquot and store the Contactin 1 antibody at -20oC.