


CRF1 Antibody (R34530)

Catalog No.	Formulation	Size
R34530-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Sheep, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1394
Applications	Immunohistochemistry (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This CRF1 antibody is available for research use only.



IHC testing of FFPE human cortex (brain) tissue with CRF1 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Corticotropin releasing hormone receptor 1, CRFR1

Application Notes

Optimal dilution of the CRF1 antibody should be determined by the researcher.

1. Reported to specifically stain mouse brain by IHC and IF.

Immunogen

Amino acids KLYYDNEKCWFGKR were used as the immunogen for this CRF1 antibody.

Storage

Aliquot and store the CRF1 antibody at -20oC.