


CRF1 Antibody (R34787)

Catalog No.	Formulation	Size
R34787-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1394
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CRF1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Corticotropin releasing hormone receptor 1; CRHR1

Application Notes

Optimal dilution of the CRF1 antibody should be determined by the researcher.

Immunogen

Amino acids NEEKKSKVHYH were used as the immunogen for this CRF1 antibody.

Storage

Aliquot and store the CRF1 antibody at -20oC.