


Cripto-1 Antibody (R34183)

Catalog No.	Formulation	Size
R34183-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55997
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Cripto-1 antibody is available for research use only.



IHC testing of FFPE mouse embryonic liver with Cripto-1 antibody. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Cryptic family 1; CFC1

Application Notes

Optimal dilution of the Cripto-1 antibody should be determined by the researcher.

Immunogen

Amino acids HFTGRYCEHDQRRSEC were used as the immunogen for this Cripto-1 antibody.

Storage

Aliquot and store the Cripto-1 antibody at -20oC.