


CRTC2 Antibody / TORC2 (R34193)

Catalog No.	Formulation	Size
R34193-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	74343
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CRTC2 antibody is available for research use only.



Western blot testing of rat spinal cord lysate with CRTC2 antibody at 0.03ug/ml.
Predicted molecular weight: ~73 kDa.

Description

Additional name(s) for this target protein: CREB regulated transcription coactivator 2

Application Notes

Optimal dilution of the CRTC2 antibody should be determined by the researcher.

Immunogen

Amino acids HGQQPYHRPIND were used as the immunogen for this CRTC2 antibody.

Storage

Aliquot and store the CRTC2 antibody at -20oC.