


CSEN Antibody (R34307)

Catalog No.	Formulation	Size
R34307-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	30818
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This CSEN antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human KCNIP3-MYC probed with CSEN antibody [1ug/ml] in Lane A and anti-MYC [1/1000] in lane C. Mock-transfected HEK293 probed with CSEN antibody [1ug/ml] in Lane B. Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: Calsenilin, Kv channel interacting protein 3; KCNIP3

Application Notes

Optimal dilution of the CSEN antibody should be determined by the researcher.

1. This CSEN antibody is specific for isoform 1.

Immunogen

Amino acids HTPLSKKEGIK were used as the immunogen for this CSEN antibody.

Storage

Aliquot and store the CSEN antibody at -20oC.