


CSF1R Antibody (R34544)

Catalog No.	Formulation	Size
R34544-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P07333
Gene ID	1436
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CSF1R antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing Human CSF1R with MYC tag with CSF1R antibody (0.5ug/ml) in Lane A and probed with anti-MYC tag (1/1000) in lane C. Mock-transfected HEK293 probed with CSF1R antibody (1ug/ml) in Lane B.

Description

Additional name(s) for this target protein: (Macrophage) Colony stimulating factor 1 receptor

Application Notes

Optimal dilution of the CSF1R antibody should be determined by the researcher.

Immunogen

Amino acids KNIHLEKKYVR were used as the immunogen for this CSF1R antibody.

Storage

Aliquot and store the CSF1R antibody at -20oC.