


## CSF3R Antibody (R35094)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35094-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                           |
|----------------------|-----------------------------------------------------------|
| Availability         | 1-3 business days                                         |
| Species Reactivity   | Human                                                     |
| Predicted Reactivity | Cow                                                       |
| Format               | Antigen affinity purified                                 |
| Clonality            | Polyclonal (goat origin)                                  |
| Isotype              | Goat Ig                                                   |
| Purity               | Antigen affinity                                          |
| Gene ID              | 1441                                                      |
| Applications         | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations          | This CSF3R antibody is available for research use only.   |



Western blot testing of peripheral blood lymphocyte lysate with CSF3R antibody at 1ug/ml. Predicted molecular weight: ~92kDa.

## Description

Additional name(s) for this target protein: Colony stimulating factor 3 receptor (granulocyte); CD114

## Application Notes

Optimal dilution of the CSF3R antibody should be determined by the researcher.

## Immunogen

Amino acids PPELGKSPLTHYT were used as the immunogen for this CSF3R antibody.

## Storage

Aliquot and store the CSF3R antibody at -20oC.