


CSK Antibody (R33934)


Catalog No.	Formulation	Size
R33934-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1445
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CSK antibody is available for research use only.



IHC staining of FFPE human colon with CSK antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.



Western blot testing of Jurkat lysate with CSK antibody at 0.5ug/ml. Predicted molecular weight: ~50 kDa.



Western blot of mouse spleen lysate with CSK antibody at 1ug/ml. Predicted molecular weight: ~50 kDa.

Description

Additional name(s) for this target protein: C-src tyrosine kinase

Application Notes

Optimal dilution of the CSK antibody should be determined by the researcher.

Immunogen

Amino acids EQLEHIKTHELHL were used as the immunogen for this CSK antibody.

Storage

Aliquot and store the CSK antibody at -20°C.