


CSNK2B Antibody (R35008)


Catalog No.	Formulation	Size
R35008-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Pig
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1460
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CSNK2B antibody is available for research use only.



Western blot testing of A) mouse NIH3T3 and human B) HeLa, C) K562, D) MCF7 and E) HepG2 lysates with CSNK2B antibody at 0.3ug/ml. Predicted molecular weight: ~25kDa.



Western blot testing of pig testis lysate with CSNK2B antibody at 0.3ug/ml. Predicted molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: Casein kinase 2, beta polypeptide

Application Notes

Optimal dilution of the CSNK2B antibody should be determined by the researcher.

Immunogen

Amino acids PEYRPKR PANQFVPR were used as the immunogen for this CSNK2B antibody.

Storage

Aliquot and store the CSNK2B antibody at -20oC.