


CX3CR1 Antibody / CMKBRL1 (R33592)

Catalog No.	Formulation	Size
R33592-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1524
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CX3CR1 antibody is available for research use only.



Western blot testing of KELLY [A] and U251-MG [B] lysates with CX3CR1 antibody at 2ug/ml. Predicted molecular weight ~40kDa.

Description

Additional name(s) for this target protein: Chemokine (C-X3-C motif) receptor 1

Application Notes

Optimal dilution of the CX3CR1 antibody should be determined by the researcher.

Immunogen

Amino acids DQFPESVTENFEYD were used as the immunogen for this CX3CR1 antibody.

Storage

Aliquot and store the CX3CR1 antibody at -20oC.