


CXCR5 Antibody | C-X-C chemokine receptor type 5 [clone 20C58] (RQ8898)

Catalog No.	Formulation	Size
RQ8898	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	20C58
Purity	Affinity chromatography
UniProt	P32302
Localization	Cytoplasm
Applications	Western Blot : 1:500 Immunohistochemistry (FFPE) : 1:50
Limitations	This CXCR5 antibody is available for research use only.



IHC staining of FFPE rat spleen tissue with CXCR5 antibody, HRP-secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) human MOLT4, 2) human HEL, 3) human Jurkat, 4) human U-2 OS, 5) rat PC-12, 6) rat C6, 7) mouse RAW264.7 and 8) mouse ANA-1 cell lysate with CXCR5 antibody. Predicted molecular weight ~43 kDa.