


CXCR6 Antibody (R35867)


Catalog No.	Formulation	Size
R35867-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10663
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.3-1ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This CXCR6 antibody is available for research use only.



Western blot testing of HEK293 lysate with CXCR6 antibody at 1ug/ml. Predicted molecular weight: ~39 kDa.



IHC testing of FFPE human thymus tissue with CXCR6 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Chemokine (C-X-C motif) receptor 6, BONZO, CD186; CXCR6

Application Notes

Optimal dilution of the CXCR6 antibody should be determined by the researcher.

Immunogen

Amino acids SEDNSKTFSASHN were used as the immunogen for this CXCR6 antibody.

Storage

Aliquot and store the CXCR6 antibody at -20oC.