


Cyb5r3 Antibody (R35083)

Catalog No.	Formulation	Size
R35083-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	109754
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Cyb5r3 antibody is available for research use only.



Western blot testing of mouse liver lysate with Cyb5r3 antibody at 0.1ug/ml. Predicted molecular weight (mouse): ~34/32kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Cytochrome b5 reductase 3, DIA1

Application Notes

Optimal dilution of the Cyb5r3 antibody should be determined by the researcher.

Immunogen

Amino acids PNLERVGHPKERC were used as the immunogen for this Cyb5r3 antibody.

Storage

Aliquot and store the Cyb5r3 antibody at -20oC.