


Cyclin B3 Antibody (CCNB3) (F52195)

Catalog No.	Formulation	Size
F52195-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52195-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8WWL7
Applications	Western Blot : 1:1000
Limitations	This Cyclin B3 antibody is available for research use only.



Description

CCNB3 is a mitotic cyclin that shares homology with cyclins of the A and B types. [OMIM]

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the Cyclin B3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1311-1338 from the human protein was used as the immunogen for this Cyclin B3 antibody.

Storage

Store at 4oC for up to one month. For long term, aliquot the Cyclin B3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.