


CYFIP1 Antibody / Cytoplasmic FMR1-interacting protein 1 [clone 29C25] (RQ8784)

Catalog No.	Formulation	Size
RQ8784	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	29C25
Purity	Affinity chromatography
UniProt	Q7L576
Applications	Western Blot : 1:500-1:2000
Limitations	This CYFIP1 antibody is available for research use only.



Western blot testing of 1) human HeLa, 2) human SH-SY5Y, 3) human U-2 OS, 4) rat brain, 5) rat C6, 7) mouse brain and 8) mouse Neuro-2a cell lysate with CYFIP1 antibody. Predicted molecular weight ~145 kDa but multiple smaller isoforms.

Application Notes

Optimal dilution of the CYFIP1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human CYFIP1 protein was used as the immunogen for the CYFIP1 Antibody.

Storage

Store the CYFIP1 antibody at -20°C.

