


CYP26A1 Antibody (R34861)

Catalog No.	Formulation	Size
R34861-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1592
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml Immunocytochemistry : 2-5ug/ml. ELISA (peptide) LOD : 1:16000
Limitations	This CYP26A1 antibody is available for research use only.



IHC testing of PFA-fixed and TX100-treated HeLa cells with CYP26A1 antibody at 2.5ug/ml (DAPI nuclear counterstaining in blue).

Description

Additional name(s) for this target protein: Cytochrome P450, family 26, subfamily A, polypeptide 1

Application Notes

Optimal dilution of the CYP26A1 antibody should be determined by the researcher.

Immunogen

Amino acids NLPARFTHFHGE were used as the immunogen for this CYP26A1 antibody.

Storage

Aliquot and store the CYP26A1 antibody at -20oC.