


## CYP7B1 Antibody (R36390)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36390-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                                          |
|--------------------|--------------------------------------------------------------------------|
| Availability       | 1-3 business days                                                        |
| Species Reactivity | Human                                                                    |
| Format             | Antigen affinity purified                                                |
| Clonality          | Polyclonal (goat origin)                                                 |
| Isotype            | Goat Ig                                                                  |
| Purity             | Antigen affinity                                                         |
| Gene ID            | 9420                                                                     |
| Applications       | Immunohistochemistry (FFPE) : 3-6ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This CYP7B1 antibody is available for research use only.                 |



IHC testing of FFPE human breast with CYP7B1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Cytochrome P450, family 7, subfamily B, polypeptide 1

### Application Notes

Optimal dilution of the CYP7B1 antibody should be determined by the researcher.

### Immunogen

Amino acids YPDSVDLFRYKVKS were used as the immunogen for this CYP7B1 antibody.

## Storage

Aliquot and store the CYP7B1 antibody at -20oC.