


## Cytokeratin 19 Antibody (R34300)

Catalog No.	Formulation	Size
R34300-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	3880
<b>Applications</b>	Western Blot : 0.03-0.1ug/ml Enzyme ImmunoAssay : used as capture/detect for recombinant protein ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This Cytokeratin 19 antibody is available for research use only.



Western blot testing of human colon cancer lysate with Cytokeratin 19 antibody at 0.03ug/ml. Predicted/observed molecular weight ~43/40-43kDa.

### Description

Additional name(s) for this target protein: Keratin 19, CK19; KRT19

### Application Notes

Optimal dilution of the Cytokeratin 19 antibody should be determined by the researcher.

### Immunogen

Amino acids HTEQLQMSRSEVTD were used as the immunogen for this Cytokeratin 19 antibody.

## Storage

Aliquot and store the Cytokeratin 19 antibody at -20oC.