


Cytokeratin 20 Antibody Biotin Conjugate (R34023BTN)

Catalog No.	Formulation	Size
R34023BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54474
Applications	Western Blot : 1-3ug/ml Enzyme ImmunoAssay : suitable as a detect ELISA (peptide) LOD : 1:128000
Limitations	This Cytokeratin 20 antibody is available for research use only.



Western blot testing of human colon lysate with biotinylated Cytokeratin 20 antibody at 1ug/ml. Predicted molecular weight: ~49kDa.

Description

Additional name(s) for this target protein: Keratin 20, CK20; KRT20

Application Notes

Optimal dilution of the Cytokeratin 20 antibody should be determined by the researcher.

Immunogen

Amino acids ETNAPRAGRDIS were used as the immunogen for this Cytokeratin 20 antibody.

Storage

Aliquot and store the Cytokeratin 20 antibody at -20oC.