


Dars Antibody (R36445)

Catalog No.	Formulation	Size
R36445-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q922B2
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) : 1:128000
Limitations	This Dars antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with Dars antibody at 1ug/ml. Predicted molecular weight ~57 kDa.

Description

Additional name(s) for this target protein: Aspartyl-tRNA synthetase

Application Notes

The stated application concentrations are suggested starting points. Titration of the Dars antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Amino acids SRKGQEKPREIVD were used as the immunogen for this Dars antibody.

Storage

Aliquot and store the Dars antibody at -20oC.