


DAX-1 Antibody / NR0B1 (RQ8276)


Catalog No.	Formulation	Size
RQ8276	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P51843
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunofluorescence : 5ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This DAX-1 antibody is available for research use only.



Immunofluorescent staining of FFPE human A549 cells with DAX-1 antibody (red) and Beta Tubulin mAb (green). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of human 1) MCF7, 2) HeLa, 3) Jurkat and 4) U-87 MG cell lysate with DAX-1 antibody. Predicted molecular weight: 44-52 kDa (two isoforms).

Description

DAX1 (dosage-sensitive sex reversal, adrenal hypoplasia critical region, on chromosome X, gene 1) is a nuclear receptor protein that in humans is encoded by the NR0B1 gene (nuclear receptor subfamily 0, group B, member 1). This gene encodes a protein that contains a DNA-binding domain. The encoded protein acts as a dominant-negative regulator of transcription which is mediated by the retinoic acid receptor. This protein also functions as an anti-testis gene by acting antagonistically to Sry. Mutations in this gene result in both X-linked congenital adrenal hypoplasia and hypogonadotropic hypogonadism.

Application Notes

Optimal dilution of the DAX-1 antibody should be determined by the researcher.

Immunogen

An E.coli-derived human recombinant protein (V126-Q404) was used as the immunogen for the DAX-1 antibody.

Storage

After reconstitution, the DAX-1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.