


Daxx Antibody [clone PCRP-DAXX-8C2] (V8935)


Catalog No.	Formulation	Size
V8935-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V8935-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V8935SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG
Clone Name	PCRP-DAXX-8C2
Purity	Protein A/G affinity
UniProt	Q9UER7
Localization	Nucleus, Cytoplasm
Applications	Flow Cytometry : 1-2ug/million cells Immunofluorescence : 1-2ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This Daxx antibody is available for research use only.




Immunofluorescent staining of PFA-fixed human HeLa cells using Daxx antibody (green, clone PCRP-DAXX-8C2).



FACS staining of PFA-fixed human HeLa cells with Daxx antibody (blue, clone PCRP-DAXX-8C2), and unstained cells (red).



IHC staining of FFPE human kidney tissue with Daxx antibody (clone PCRP-DAXX-8C2). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free Daxx antibody (clone PCRP-DAXX-8C2) as confirmation of integrity and purity.

Description

Daxx, or Death domain-associated protein 6, is a transcription corepressor known to repress transcriptional potential of several sumoylated transcription factors. [UniProt]

Application Notes

Optimal dilution of the Daxx antibody should be determined by the researcher.

Immunogen

Recombinant full-length human Death domain-associated protein 6 protein was used as the immunogen for the Daxx antibody.

Storage

Aliquot the Daxx antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.