


DAZL Antibody (R34115)

Catalog No.	Formulation	Size
R34115-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1618
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DAZL antibody is available for research use only.



Western blot testing of rat testis lysate with DAZL antibody at 1ug/ml. Predicted molecular weight: ~32 kDa.

Description

Additional name(s) for this target protein: Deleted in azoospermia-like

Application Notes

Optimal dilution of the DAZL antibody should be determined by the researcher.

Immunogen

Amino acids GNGPQKKSVDR were used as the immunogen for this DAZL antibody.

Storage

Aliquot and store the DAZL antibody at -20oC.