


DBC1 Antibody (R33992)


Catalog No.	Formulation	Size
R33992-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57805
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This DBC1 antibody is available for research use only.



Western blot testing of HEK293 [A] and MCF7 [B] nuclear lysates with DBC1 antibody at 0.5ug/ml. Predicted molecular weight: ~103/140kDa (unmodified/acetylated).



Western blot testing of mouse NIH3T3 lysate with DBC1 antibody at 2ug/ml. Predicted molecular weight: ~103/140kDa (unmodified/acetylated).

Description

Additional name(s) for this target protein: CCAR2, Cell cycle and apoptosis regulator 2 ; CCAR2

Application Notes

Optimal dilution of the DBC1 antibody should be determined by the researcher.

Immunogen

Amino acids ESHNRFSATEVTNK were used as the immunogen for this DBC1 antibody.

Storage

Aliquot and store the DBC1 antibody at -20oC.