


DBNL Antibody (R33109)


Catalog No.	Formulation	Size
R33109-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	28988
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DBNL antibody is available for research use only.



Western blot testing of Jurkat lysate with DBNL antibody at 1ug/ml. Predicted molecular weight: ~49kDa.



Western blot testing of NIH3T3 lysate with DBNL antibody at 1ug/ml. Predicted molecular weight: ~49kDa.

Description

Additional name(s) for this target protein: Drebrin-like; HIP-55

Application Notes

Optimal dilution of the DBNL antibody should be determined by the researcher.

Immunogen

Amino acids AANLSRNGPALQE were used as the immunogen for this DBNL antibody.

Storage

Aliquot and store the DBNL antibody at -20oC.