


Dcdc2a Antibody (R35120)

Catalog No.	Formulation	Size
R35120-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51473
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Dcdc2a antibody is available for research use only.



Western blot testing of mouse brain lysate with A) Dcdc2a antibody, cat # R35120 (0.1ug/ml) and B) [DCDC2 antibody, cat # R33450](#) (0.1ug/ml). Predicted molecular weight (mouse): ~52/41kDa (isoforms 1/2), observed here at ~30kDa.

Description

Additional name(s) for this target protein: Doublecortin domain containing 2a

Application Notes

Optimal dilution of the Dcdc2a antibody should be determined by the researcher.

Immunogen

Amino acids PQPLKRKGKKDSN were used as the immunogen for this Dcdc2a antibody.

Storage

Aliquot and store the Dcdc2a antibody at -20oC.