


DCN1 Antibody (R35717)


Catalog No.	Formulation	Size
R35717-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54165
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DCN1 antibody is available for research use only.



Western blot testing of HepG2 lysate with DCN1 antibody at 0.5ug/ml. Predicted molecular weight: ~30kDa.



IHC staining of FFPE human kidney with DCN1 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Defective in cullin neddylation 1, domain containing 1; DCUN1D1

Application Notes

Optimal dilution of the DCN1 antibody should be determined by the researcher.

Immunogen

Amino acids RPQIAGTKSTT were used as the immunogen for this DCN1 antibody.

Storage

Aliquot and store the DCN1 antibody at -20oC.