


DCXR Antibody (R34140)


Catalog No.	Formulation	Size
R34140-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51181
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DCXR antibody is available for research use only.



Western blot testing of human kidney lysate with DCXR antibody at 0.1ug/ml. Predicted molecular weight: ~26kDa.



IHC staining of FFPE human kidney with DCXR antibody at 1ug/ml. HIER: microwaved with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Dicarbonyl Reductase

Application Notes

Optimal dilution of the DCXR antibody should be determined by the researcher.

Immunogen

Amino acids GSTLPVEGGFWAC were used as the immunogen for this DCXR antibody.

Storage

Aliquot and store the DCXR antibody at -20oC.