


DEF6 Antibody (R33147)

Catalog No.	Formulation	Size
R33147-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	50619
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This DEF6 antibody is available for research use only.



Western blot testing of human thymus lysate with DEF6 antibody at 1ug/ml. Predicted molecular weight: ~74kDa.

Description

Additional name(s) for this target protein: Differentially expressed in FDCP 6

Application Notes

Optimal dilution of the DEF6 antibody should be determined by the researcher.

Immunogen

Amino acids ALRKELLKSIWYA were used as the immunogen for this DEF6 antibody.

Storage

Aliquot and store the DEF6 antibody at -20oC.