


Desmoglein 1 Antibody (R36403)

Catalog No.	Formulation	Size
R36403-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q02413
Gene ID	1828
Localization	Cell surface, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Desmoglein 1 antibody is available for research use only.



IHC testing of FFPE human skin with Desmoglein antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: DSG1; DSG; DG

Application Notes

Optimal dilution of the Desmoglein 1 antibody should be determined by the researcher.

Immunogen

Amino acids YTTSDLTKPSVHVHD were used as the immunogen for this Desmoglein 1 antibody.

Storage

Aliquot and store the Desmoglein 1 antibody at -20oC.