


DLL4 Antibody (Center Region) (R34555)

Catalog No.	Formulation	Size
R34555-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54567
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This DLL4 antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] lung lysates with DLL4 antibody at 0.3ug/ml. Predicted/observed molecular weight 75–85kDa.

Description

Additional name(s) for this target protein: Delta ligand 4

Application Notes

Optimal dilution of the DLL4 antibody should be determined by the researcher.

Immunogen

Amino acids KNTNQKKELEVDC were used as the immunogen for this DLL4 antibody.

Storage

Aliquot and store the DLL4 antibody at -20oC.