


DMP1 Antibody (R33900)

Catalog No.	Formulation	Size
R33900-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1758
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This DMP1 antibody is available for research use only.



Western blot testing of human kidney lysate with DMP1 antibody at 1ug/ml. Predicted molecular weight: ~56kDa.

Description

Additional name(s) for this target protein: Dentin matrix acidic phosphoprotein

Application Notes

Optimal dilution of the DMP1 antibody should be determined by the researcher.

Immunogen

Amino acids ENSEQANTQDSG were used as the immunogen for this DMP1 antibody.

Storage

Aliquot and store the DMP1 antibody at -20oC.